Agilent 3200P pH Meter pH 计

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Agilent Technologies

Notices

© Agilent Technologies, Inc. 2012

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

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Warranty

The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

U.S. Government Restricted Rights. Software and technical data rights granted to the federal government include only those rights customarily provided to end user customers. Agilent provides this customary commercial license in Software and technical data pursuant to FAR 12.211 (Technical Data) and 12.212 (Computer Software) and, for the Department of Defense, DFARS 252.227-7015 (Technical Data -Commercial Items) and DFARS 227.7202-3 (Rights in Commercial Computer Software or Computer Software Documentation).

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

Contents

1	Installation	5
	Tools and Components Needed for Installation	6
	Installing the 3200P pH Meter	7
	Installing the Power Adaptor	9
	Installing the Ground Line	10
	Installing the Electrochemical Data Collecting Software	10
	More Information	10
2	Operation	11
	Introduction	12
	Safety Information	14
	Agilent Customer Portal	15
	Term Definitions	16
	Features	17
	Specifications	18
	The Operating Panel	19
	Display	19
	Keypad	19
	Instrument control	21
	Setting Up the pH Meter	22
	To start up the 3200P pH Meter	22
	Accessing setup parameters	23
	To manually set the solution temperature	27
	To set up the system (System Setup)	29
	To select the display language	31
	To set the meter to turn off automatically (Auto Power Off)	31

Using the pH Meter	32
To prepare an electrode for use	32
To adjust mV zero	33
To calibrate pH electrode (Calib pH EC)	36
To measure pH value	39
To measure mV value	40
To measure temperature	40
To turn off the 3200P pH Meter	40

Storing Data	41
To delete saved data	41
To output data	41

Viewing Parameters During Measurement	42
To view saved data	43
To view last calibrated data	45

Software Operation	47
To install the software	47
To run software	47


3 Troubleshooting and Maintenance 49

Troubleshooting Common Problems	50
---------------------------------	----

General Troubleshooting Procedure	52
Check the temperature measurement	52
Check the potential (mV) value	53
Check the pH value	53
Restore factory default settings	54

Meter Self-Diagnostics and Messages	56
Severe errors	56
Error warning	57

Maintenance	59
Electrode cleaning	59
Electrode storage	59

1

Installation

Tools and Components Needed for Installation	6
Installing the 3200P pH Meter	7
Installing the Power Adaptor	9
Installing the Ground Line	10
Installing the Electrochemical Data Collecting Software	10
More Information	10

Use the tools and accessories shipped with 3200P pH Meter during installation.

Tools and Components Needed for Installation

Agilent provides the special tools needed for installation. You can find the following parts in the shipping case of the pH Meter (Meter Kit).

- Electrode holder (G4389A)
- Electrode, for example P3211 pH combination electrode (5190-3988)
- Power adapter (5185-8389)

Installing the 3200P pH Meter

Open the meter kit package case. Remove the pH Meter, multifunction electrode holder, and other parts.


Figure 1 Installing the Electrode holder

The sockets for the pH electrode and ATC probe are on the back panel of the meter. Insert the connectors of the pH electrode and ATC probe (5190-3998) into their respective sockets. See [Figure 6](#) on page 13. Run the electrode cables along the electrode holder as shown in [Figure 2](#) on page 8.

1 Installation

Figure 2 Installing the electrodes

Installing the Power Adaptor

The universal power adapter that is included with the meter is the only power adapter recommended for use with this unit.

The external electrical power adapter is rated at 100 to 240 VAC, 1 A, 50/60 Hz.

This power adaptor provides a plug converter which supports several kinds of plugs. Choose the appropriate power plug. Connect the output plug of the power adapter to the power input on the meter. A click will be heard when the plug is properly engaged.

Figure 3 Installing the Power Adaptor

Installing the Ground Line

The meter ships with an optional ground line. During use, if another device (such as a constant temperature bath) causes electrical interference and an unstable reading, install the ground line. Connect one terminal of the ground line to the meter ground socket () and the other terminal to the interference source.


Installing the Electrochemical Data Collecting Software

The optional G4390A software provides communications between the meter and a computer. To input data to a computer, install the Electrochemical Data Collecting Software on a computer. Connect the computer to the meter with a USB cable. See the G4390A software documentation for more information.

If not using the G4390A software, Agilent provides downloadable data printing software on the Agilent Customer Portal (see “[Agilent Customer Portal](#)” on page 15). To use this software, download it from the portal and install it. Then connect the meter to the PC using the USB cable. Refer to the data printing software documentation for more information

More Information

For more information, visit the Agilent website at <http://www.agilent.com/chem>. See also the instructions provided with each electrode.

2 Operation

Introduction	12
Safety Information	14
Agilent Customer Portal	15
Term Definitions	16
Features	17
Specifications	18
The Operating Panel	19
Setting Up the pH Meter	22
Using the pH Meter	32
Storing Data	41
Viewing Parameters During Measurement	42
Software Operation	47

This chapter describes how to use the 3200P pH Meter.

Introduction

The 3200M pH Meter is suitable for accurately measuring the pH value of various solutions.

During measurement, the LCD display shows the current state of the 3200P pH Meter. During measurement, you can also change the displayed parameters.

The 3200P pH Meter is composed of a meter body and an electrode system. The electrode system shown in [Figure 4](#) includes a pH combination electrode and an ATC probe.

See [Figure 4](#), [Figure 5](#), and [Figure 6](#).

Figure 4 Major parts of the pH 3200P system

Figure 5 Front view of the pH 3200P Meter


Figure 6 Back view of the pH 3200P Meter

Safety Information

Before operating your meter, be sure to read the installation and operation instructions.

The most common safety issues when working with the meter include:

- Using any power source other than the original power adaptor may produce some safety problems.
- Be sure that the meter has a good ground connection.
- Prevent exposure to corrosive gases.
- Keep the sockets of the meter clean and dry. Avoid contact with acidic, alkaline, and salt solutions.
- The meter can be used continuously for a long time. After a measurement, soak all electrodes in distilled water.
- To protect high resistance components from damage, install the short circuit plug (G4383-40000) onto the pH electrode socket when the meter is not connected to a pH electrode. (The short-circuit plug ships installed in the pH socket of the meter, see [Figure 6](#).) Keep the short-circuit plug clean and dry. The plug can be damaged by humidity. A damaged short-circuit plug can damage the meter.

Agilent Customer Portal

Agilent also provides customized information for the products you own through a customer portal. This web service provides many customizable services as well as information related directly to your Agilent products and orders. Log onto the portal at ***<http://www.agilent.com/chem>***.

Term Definitions

pH Slope	The mV variation of each unit of pH represented by mV/pH or percent (%) theoretical slope
pH E_0	Also called <i>Zero Potential</i> and is the mV value at pH 7
One-point pH Calibration	Calibration with one pH buffer
Multi-point pH Calibration	Calibration with more than one pH buffer

Features

The 3200P pH Meter is a state-of-the-art and customer-friendly benchtop analytical instrument, designed for the measurement of pH, potential (mV), and temperature in aqueous solutions. It features:

- Measurement of pH, potential (mV), and temperature
- Automatic pH buffer recognition for buffers prepared according to NIST, DIN and GB standards
- Supports multi-point calibration (up to five point calibration)
- LCD display
- Supports GLP:
 - Requires setting an operator number
 - Records all measurements
 - Records, views and prints calibration data
 - All saved data meets GLP. Stores a maximum of 200 sets of measurement data
- Can view, print and delete stored measurement data
- Supports three measuring modes: Continuous Mode, Timed Reading, and Auto-Lock Mode
- A USB port can be connected to a PC when using the optional G3490A software.
- When the meter is powered off, the stored measuring data, calibrated data, and setting parameters will not be lost.
- Can turn itself off to conserve power
- Provides self-diagnostic messages and functions to support troubleshooting issues that occur
- Provides LCD back lighting for use under dim lighting conditions
- Meter firmware can be updated for any new features or bug fixes that arise.

Specifications

The major specifications of the 3200P pH Meter include measuring range, resolution, accuracy of the electronics, normal working conditions, size and weight.

Measurement range

pH: -2.000 to 20.000
mV: -1999.9 mV to 1999.9 mV
Temperature: -5.0 to 110.0 °C

Resolution

pH: $\pm 0.1/0.01/0.001$ pH/pX
mV: ± 0.1 mV
Temperature: ± 0.1 °C

Accuracy

pH: ± 0.002
mV: $\pm 0.03\%$ full scale
Temperature: ± 0.1 °C

Normal operating conditions

Ambient (environment) temperature: 0 to 40 °C
Corresponding humidity: $\leq 85\%$
Power supply: Power adaptor (5185-8389)
Power input: 100 to 240 VAC, 1 A
Power output: 9 VDC, 1 A
No performance-affecting vibrations nearby
No corrosive gases in ambient air
No strong electromagnetic interference
Size (mm): 190 length \times 190 width \times 105 height
Weight (kg): approximately 1 kg

The Operating Panel

Display


Figure 7 The display

The display shows the working condition and current setting of the 3200P pH Meter. After the power supply is connected, press **[On/Off]** to turn on the meter. After a successful power up, the meter model, version number, and other information displays. After the system self-check, the meter goes to the initial state as shown in [Figure 7](#). The left side of the screen shows current system time and the right side shows the current measuring mode, parameters, and the last calibration data.

Keypad

The meter has 15 operating keys. The keys are shown in [Figure 8](#) on page 20 and described in [Table 1](#) on page 20. Most keys have two values, a numeric entry and a function. Which one applies depends on the display; if it is waiting for a number, the numeric entry applies, otherwise, the function applies. To complete a numeric entry, press **[Enter]**.


Figure 8 Operating keys

Table 1 Key values

Key symbol	Numeric entry	Function
1/Output	1	Output presently displayed data
2/▲	2	Move up or select mode
3/Save	3	Store measuring data
4/◀	4	Move left
5/Setup	5	Set different function
6/▶	6	Move right
7/View	7	View stored or calibrated data
8/▼	8	Move down or select mode
9/Mode	9	Switch display window or parameter in measuring state
0/Measure	0	Begin measurement in initial state
–/Delete	– (negative number)	Delete function Can delete saved data when you view them.
./Calibrate	. (decimal point)	Calibrate electrode
Enter		Enter
Cancel		Cancel
On/Off		Turn meter on or off

Instrument control

The 3200P pH Meter is usually controlled by the keypad. It can also be controlled through a computer using the optional G4390A software. For details on how to use the optional data collecting software, see the information provided with the software.

Setting Up the pH Meter

This section describes the tasks that an operator will perform when starting up the 3200P pH Meter.

Using the meter involves:

- 1 Selecting the measurement mode and the display parameters.
- 2 Preparing the electrode.
- 3 Measuring the pH.

To start up the 3200P pH Meter

Always operate the meter in an environment that meets the requirements listed in “[Normal operating conditions](#)” on page 18.

- 1 Connect the meter to the power supply.
- 2 Ensure the meter has a good ground.
- 3 Check the electrode connector and meter socket. The electrode must have a reliable connection.
- 4 Press **[On/Off]** to turn the meter on.

The initial state of the meter is shown in [Figure 9](#). The left side shows the system time. The right side shows the setting, measuring mode, pH parameter, and the last calibration data.


Figure 9 Function setup in the initial state

Accessing setup parameters

From the initial state display, press **[Setup]** to display a menu with choices to define system parameters (date, time, temperature, and so forth) or select the measuring mode. Choices include: **Measuring Mode**, **System Setup**, **Language select**, **Set Manual Temp**, and **Set Condition**. See [Figure 10](#).


Figure 10 The setup menu

To make a selection, use the arrow keys to highlight the desired parameter, then press **[Enter]** to view its settings. Press **[Cancel]** to exit.

[Measuring Mode]

Select to set the measuring mode used when taking measurements: Continuous Mode, Timed Reading (Timed Reading Mode), Auto-Lock Mode (Auto-Lock Reading Mode). See [“To select a measuring mode”](#) on page 25.

[Set Manual Temp]

Select to manually enter the current sample solution temperature. If the meter is not connected with an ATC probe, the temperature value must be input manually. Otherwise, the meter uses the ATC probe value. See [“To manually set the solution temperature”](#) on page 27.

[Set Condition]

Select this mode to enter an Auto-Lock condition in Auto-Lock reading mode. See “[To set a condition for Auto-Lock mode](#)” on page 28.

[System Setup]

Set up system time and GLP Standard. See “[To set up the system \(System Setup\)](#)” on page 29.

[Language Select]

Supports Chinese and English version. See “[To select the display language](#)” on page 31.

[Auto Power Off]

Supports the function of turning off automatically. See “[To set the meter to turn off automatically \(Auto Power Off\)](#)” on page 31.

[Set Default]

Select to delete all custom settings and restore the factory default settings. See “[To restore factory defaults \(Set Default\)](#)” on page 54.

To select a measuring mode

The meter supports three measuring modes: **Continuous Mode**, **Timed Reading**, and **Auto-Lock Mode**. After selecting a parameter and measuring mode, the next measurement will use the new settings.


Figure 11 Measuring mode

To select a measuring mode:

- 1 Press **[Setup]** and then **[Enter]**. The meter displays is similar to [Figure 11](#). The left side displays parameters. The right side lists the measuring modes. A ✓ next to a mode means that this mode has been selected in the current setup.
- 2 Use the arrow keys to highlight the desired mode.
- 3 Press **[Setup]** to select the highlighted mode.
- 4 Press **[Enter]** to save the new setup, exit the **Setup** mode, and return to the initial state.

Press **[Cancel]** to exit the **Setup** mode without saving changes.

Measuring modes

Continuous mode

This is the most commonly used measuring mode. The meter performs the measurement and calculates and displays the result continuously. You can view the calibrated parameter, calibrate the electrode, and save or output the result during the measurement. After measurement, press **[Cancel]** and then **[Enter]** to exit the measuring mode.

Timed Reading mode

In this mode, the meter continuously measures the pH for a specified period of time. For example, to measure pH continuously for 30 minutes, select this mode and set the time interval to **30**. When measurement begins, the meter continuously displays the results of measurements and calculations. When the set time period ends, the meter saves the measuring data. (If the USB interface is connected to a PC, the meter also prints out the measurement data). Then the meter begins the next measurement automatically. To stop measurement, press **[Cancel]** and then **[Enter]**.

In this mode, you enter the desired time interval. The default value for the time interval is 10 minutes. The range is 1 to 99 minutes.

Auto-Lock Reading mode (Auto-Lock mode)

Auto-Lock mode requires that you set a condition as described in “[To set a condition for Auto-Lock mode](#)” on page 28. When measurement begins, the meter continuously displays the results of measurements and calculations until the measurements meet the preset Auto-Lock conditions, at which point the measurement ends.

During a measurement in this mode, you can view measuring parameters, calibrate the electrode, and so forth. After measurement ends, you can save and output the results.

Press [**Cancel**] to exit this mode or press [**Measure**] to begin the next measurement.

To manually set the solution temperature

If using an ATC probe (an ATC probe is connected to the meter), the meter uses the temperature data from the probe. Otherwise, the user must enter the temperature value of the sample solution.


Figure 12 Entering temperature

- 1 Press [**Setup**] and select **Set Manual Temp**.
- 2 Enter the temperature in the dialog.

To set a condition for Auto-Lock mode

In Auto-Lock mode (only), **Set Condition** is used to set the termination conditions for measuring the pH. The meter measures the pH until the condition is met and is stable for the duration specified in **Auto-Lock Time**.

For example, [Figure 13](#) shows that the condition for pH is **0.010pH**. When the pH measurement stabilizes within ± 0.01 pH, and remains stable for the **Auto-Lock Time** of 5 seconds, measurement ends.


Figure 13 Setting a condition

The valid time range is 1 to 200 seconds.

To set up the system (System Setup)

System setup includes the following items:

Calib Interval	Calib Interval is the time interval between recalibration prompts. The meter counts the time since the last calibration and if the time counted exceeds the entered calibration time interval, the meter displays a popup window to remind the user to recalibrate the electrode. Calib Interval is measured in hours. Input zero (0) hours to disable.
Operator No	Operator No is a three-digit number in the range 000 to 200. Every set of recorded data contains the Operator No. field. This field is used to record the operator performing measurements for regulated (GLP) laboratories. If you need to track the operator, assign a unique number for each operator.
Time	System time and date. See “To set the date and time” on page 30.
Auto Delete	Enable to automatically overwrite the oldest stored data entries with new data if the data storage is full.

- 1 Press **[Setup]**.
- 2 Select **System Setup**, and press **[Enter]** to enter **System Setup**. See [Figure 14](#).


Figure 14 System setup

- 3 Use the arrow keys to select the desired parameter.
- 4 Press **[Setup]** to modify the corresponding parameter.
Enter the new setting using the keypad.
- 5 Press **[Enter]** to save the change, exit **System Setup**, and return to the initial state.

CAUTION

The Auto Delete function is designed for cases when stored data reaches the limit of the memory of the meter, which allows you to choose whether to automatically overwrite the old data with the new data. For example, the meter stores up to 200 sets of data. When the user wants to store the 201st set of data and Auto Delete is checked, the meter deletes the first set of data and stores the 201st data. If **Auto Delete** is disabled, the new data will not be stored.

If the **Auto Delete** function is disabled, the meter will remind the user to store new data when the meter runs out of memory

To set the date and time

- 1 From the system setup screen, select **Time**.
- 2 Press **[Setup]** to open a data entry window as shown in [Figure 15](#) on page 31.


Figure 15 Set date and time

This window displays current time setup including year, month, day, hour, minute, and second. To modify the time setup, use the arrow keys to highlight the proper time cell, then press **[Setup]** and input the correct value. Press **[Enter]** to save. When all changes have been made, press **[Enter]** to save the final setup. Press **[Cancel]** to exit from **System Setup**.

To select the display language

The meter can display either Chinese or English text. To select the preferred language, from the initial state press **[Setup]**. Choose **Language Select** and then select either Chinese or English.

To set the meter to turn off automatically (Auto Power Off)

This meter can automatically turn off after a specified period of time. You can set a time period for **Auto Power Off** from 10–480 minutes. When the meter has run for the length of time that equals to **Auto Power Off** setpoint, the meter will be forced to be power off.

To disable **Auto Power Off**, set the time period to zero (0). When the meter is controlled using the optional G4390A software, **Auto Power Off** is disabled.

Using the pH Meter

For best accuracy, perform the following procedures before making measurements.

To prepare an electrode for use

A new electrode should be stored wet in the reference filling solution (5190-0545) and should be in working condition when received. Used electrodes should be stored as described in “[Electrode storage](#)” on page 59.

- 1 Check the electrode for any mechanical damage. If the storage bottle is dry, soak the electrode in reference filling solution for at least 2 hours before use.
- 2 Take the storage bottle off and keep it for future use.
- 3 Check the measuring tip for solidified electrolyte. If present, rinse off with distilled or deionized water.
- 4 Unplug the filling hole and add new reference filling solution if necessary. Fill to at least 20 mm higher than the level of sample. The liquid junction must be immersed completely in solution.
 - To guarantee the movement of reference filling solution during measurement, keep the filling hole open during measurement.
- 5 Hold the electrode measuring tip downwards and swing it several times to remove air bubbles near the sensitive glass bulb.
- 6 Install the electrode into the electrode holder as shown in [Figure 2](#) on page 8. Put the electrode in the solution with the measuring tip downwards.

Do not remove the short-circuit plug until after adjusting mv zero.

To adjust mV zero

Although the meter compensates for many factors which may affect measurement (such as temperature), it still can not guarantee the ideal condition of zero drift. For accurate measurement, perform **Adjust mV Zero** before measurement.

To adjust mV zero:

- 1** Turn off the meter.
- 2** Adjusting the zero requires the short-circuit plug (G4383-40000). If needed, disconnect the electrode and install the short circuit plug into the pH socket.
- 3** Turn on the meter and let stabilize for 30 seconds.
- 4** Press **[Setup]**, select **Adjust mV Zero**, and press **[Enter]**. The meter will display **Adjust mV zero?**.
- 5** Press **[Enter]** and the meter will calibrate mV zero.

If the potential (mV) value deviates far from zero, double check the short circuit plug connection.

After adjusting the mV zero, connect the electrode to the meter. Store the short-circuit plug in a clean, dry place. Damage to the plug can damage the meter.

To select the pH Buffer Group

The meter provides auto-recognition, to recognize several buffer groups which are prepared according to either the NIST standard, the DIN standard, or the GB standard. Each group contains several buffers.

Table 2 Buffer groups

Buffer pH	NIST group	DIN group	GB group
1.677	●		
1.680		●	●
3.557		●	
3.559			●
3.775		●	
4.003			●
4.008	●	●	
6.864	●		●
6.865		●	
7.000	●	●	
7.409			●
7.416	●	●	
9.182			●
9.184		●	
10.014	●	●	
12.454		●	
12.460			●
12.469	●		

For each Buffer Group Setup, you may select up to five buffers. Since the pH range of a buffer may overlap with the pH range of another buffer in the same buffer group, the meter limits the selection of two or more buffers with overlapping pH ranges.

To select a buffer group:

- 1 View the last pH calibrated data. See [“To view last calibrated data”](#) on page 45.
- 2 Press **[Setup]** to enter **Group Setup** as shown in [Figure 16](#).

Figure 16 Group setup

[Figure 16](#) indicates that the current Buffer Group Setup meets the DIN standard. The upper area displays the three buffer groups and the bottom area displays the buffers included in the selected buffer group.

- 3 Use the arrow keys to highlight the target buffer group and press **[Setup]**. The display is shown in [Figure 17](#).


Figure 17 Modifying a buffer

It shows icons corresponding to buffers in the NIST Buffer Group. Icons with √ indicate a selected buffer. Icons without √ indicate that this buffer is available but not selected. The bottom part of the display shows the pH value and the range of the highlighted buffer.

- 4 Use the arrow keys to highlight the target buffer, then press **[Setup]** to select or deselect it.

For example, to select the 4.008 pH buffer, use the arrow keys to highlight the icon with 4.008 pH. Press **[Setup]** and this icon is selected.

To prevent selecting buffers with overlapping pH ranges, select only buffers necessary for measurement and clear selection of other buffers.

The meter allows selection of up to five buffers for each Buffer Group Setup.

To calibrate pH electrode (Calib pH EC)

Before each measurement, Agilent recommends recalibrating the electrode. Calibration overwrites the previous calibration data.

Recalibrate the electrode regularly. To set a calibration reminder, see “[To set up the system \(System Setup\)](#)” on page 29.

To perform the calibration:

- 1 If not connected, connect the electrode to the meter. Store the short-circuit plug in a clean, dry place.
- 2 View that the desired buffers have been set up in the desired buffer group.
- 3 The display for pH calibration is shown in [Figure 18](#) on page 37. The upper area shows the current pH, potential and temperature value (in this example the slope is set to 100.00%). The bottom area shows the current calibration data (if any).


Figure 18 Calibration

- 4 Prepare one to five buffers. They can be purchased or can be prepared by you. Place them in a constant temperature environment to equilibrate.
- 5 View the pH calibration data. See [“To view last calibrated data”](#) on page 45.
- 6 Press **[Setup]** and select the pH buffer group. See [“To select the pH Buffer Group”](#) on page 34.
- 7 Return to the **Setup** window. Scroll to **Set Config Type** and press **[Enter]**.


Figure 19 Calibration setup

- 8 If using standard buffers from within a group, select **Auto-recognition**. Matched with the proper Select Group, the meter will automatically recognize these buffers. If the

meter can not recognize a buffer, it will prompt with **Calibration error**.

If you have created and are using a custom buffer, select **Manual Recognition**. You will manually input the pH of the selected buffer at the selected temperature.

- 9 Return to the **Setup** window.
- 10 Clean the pH electrode, reference electrode, and ATC probe. Place them in the buffer to be calibrated.
- 11 Press [**Measure**]. The meter will be ready to take a calibration reading for the first buffer solution selected in the buffer group.
- 12 When the reading becomes stable, press [**Enter**]. The meter will display **Storing.....** and save the calibration data.
- 13 After a few seconds, the meter will display a **Continue?** prompt window. To calibrate other buffers, select the next buffer to continue. Then repeat the measurement and save procedures. At any time during calibration, press [**Cancel**] to end calibration.

For example, to calibrate the electrode with 4.003 pH and 9.182 pH buffers, you must have selected 4.003 pH and 9.182 pH buffers in the Buffer Group Setup, otherwise the meter will not automatically recognize these two buffers. Even worse, it may lead to the wrong calibration.

The meter supports up to five calibration points. When the fifth buffer is calibrated, the meter will prompt to end calibration.

For buffers such as 6.864 pH and 7.000 pH which have overlapping pH ranges, we recommend the following calibration method.

- When you calibrate the 6.864 pH buffer, set the buffer group to contain only 6.864 pH, then calibrate it. Reset the buffer group to contain only 7.000 pH, then calibrate it.
- You could also use the manual recognition method, in which you manually input the pH of the selected buffer at the selected temperature.

You can use both auto- and manual recognition to calibrate the electrode to the buffer solutions in use.

To measure pH value

NOTE

If sure that the meter has been recently calibrated and zeroed, you can simply start taking measurements.

- 1 Before measuring pH, adjust mV zero. See “To adjust mV zero” on page 33.
- 2 Calibrate the meter. See “To calibrate pH electrode (Calib pH EC)” on page 36.
- 3 From the initial state display, press **[Setup]** and select a measuring mode. See “To select a measuring mode” on page 25.
- 4 Press **[Measure]** to start taking a measurement according to the selected mode. See Figure 20. The upper area of the screen shows the current measuring mode and system time. The central area of the screen shows pH, potential and temperature.


Figure 20 Enter continuous mode

During measurement, you can recalibrate the electrode, set parameters, and select the pH display resolution. After measurement, press **[Save]** to save data, **[Output]** to output data, or **[Cancel]** to end measurement.

To measure mV value

In any measuring mode, the meter always displays the current pH and mV value.

To measure temperature

In any measuring mode, the meter will display the current temperature value.

To turn off the 3200P pH Meter

After using the meter, ensure you have pressed [**Save**] to save data. Then, press [**On/Off**] to turn off the meter.

When not in use, all electrodes should be soaked in distilled water. If the meter is out of use for a long time, do the following:

- Disconnect the power adaptor from the power line to avoid damaging the power adaptor and the meter.
- The socket of the meter must be kept clean and dry. Avoid contact of the socket with acidic, alkaline, and salt solution.
- If disconnecting the pH electrode, install the short-circuit plug onto the meter to prevent damage to the meter.
- Store the electrode in reference filling solution. For refillable electrodes, plug the filling hole. If the electrode is out of use for a long time, put it back in the box and store it at ambient temperature.

Storing Data

The meter can save up to 200 sets of pH data. When the meter saves data depends on the measuring mode. In **Continuous Mode** and **Auto-Lock Mode**, press **[Save]** to save data when the readings become stable. In **Timed Reading** mode, the meter will automatically save measurement data after the measurement time expires. (Data can always be saved manually.)

To delete saved data

The meter can both save and delete data. Inaccurate or unwanted data can be deleted individually or all at once.

- 1 From the initial or measuring mode, press **[View]** to display data.
- 2 Scroll to the data to delete, then press **[Delete]**.

To output data

There are two ways to output current measurements, calibration data or saved data. One method is to use optional G4390A the data collecting software with the meter.

The other method is to use the free Data Printing Software.

- 1 Connect the meter to the PC with a USB cable, and run the software.
- 2 Press the **[Output]** key during measurement, while viewing calibration or stored data. This allows the user to print the corresponding data. Refer to the software operating instructions for details.

See “[Installing the Electrochemical Data Collecting Software](#)” on page 10.

Viewing Parameters During Measurement

During a measurement, the selected parameter automatically displays in the measurement window.

Figure 21 shows the measuring mode on the screen displayed with the parameter and measuring mode the user has chosen.


Figure 21 Continuous mode display

You can also view other parameters during a measurement. To view other parameter values, press **[Mode]** and the display will highlight the measuring window (see Figure 22). Press **[4/◀]** or **[6/▶]** repeatedly to switch parameters. When you stop pressing **[4/◀]** or **[6/▶]** for more than a few seconds, the meter automatically returns to the normal measuring mode display.


Figure 22 Switching parameters

To view saved data

The 3200P pH Meter lets you view the last calibrated data and stored data.

- 1 From the initial state display or measuring mode, press **[View]** and select **View Saved pH** as shown in Figure 23.

Figure 23 View menu

- 2 For saved data, use the scroll keys to select the data to view. Figure 24 shows example saved data:


Figure 24 View saved data menu

The upper area displays the current viewing mode and actual stored amount. Each page can display 10 saved data. The display format varies with different viewing modes. The data displayed include **Save time** and **Operator No.**

- 3 If desired, press **[Output]** or **[Delete]** to print or delete data respectively. The output format is shown in Figure 25 on page 44. To access the output for either of these functions, the meter must first be connected to the PC.

```

=====
MODEL                      3200P
                           PH METER
SOFT VERSION              VER 1.00
PRINT TIME                10:25:42
                           2010/06/20
OPERATOR NO               000
*****
STORED NUM:              003
*****
                           NO:001
OPERATOR NO:             000
STORED TIME:             10:19:00
                           2010/06/20
SLOPE:                   99.99%
E0:                      -0.0mV

POTENTIAL:              -68.5mV
pH:                     8.159pH
TEMP:                   25.0c
TC.TYPE:                MTC
*****
                           NO:002
OPERATOR NO:             000
STORED TIME:             10:19:00
                           2010/06/20
SLOPE:                   99.99%
E0:                      -0.0mV

POTENTIAL:              -68.5mV
pH:                     8.159pH
TEMP:                   25.0c
TC.TYPE:                MTC
*****
                           NO:003
OPERATOR NO:             000
STORED TIME:             10:19:00
                           2010/06/20
SLOPE:                   99.99%
E0:                      -0.0mV

POTENTIAL:              -68.5mV
pH:                     8.159pH
TEMP:                   25.0c
TC.TYPE:                MTC
=====

```

Figure 25 Output format

To view last calibrated data

From the initial state display or measuring mode, press **[View]** and select **pH Last Calib** (last calibrated data) as shown in Figure 26.


Figure 26 View pH Last Calib Data

Figure 27 shows an example.


Figure 27 View pH Last Calib Data

The upper area shows the last calibrated data. The bottom area shows the current buffer groups. Press the arrow keys to view detailed calibrated data.

To print current calibration data, connect the meter to the PC with a USB cable and press **[Output]**. See also “To output data” on page 41. The printout format is shown in Figure 28 on page 46.

```
=====
MODEL                                3200P
                                PH METER
VERSION                             VER 1.00
PRINT TIME                          10:11:27
                                2010/06/20
OPERATOR NO                         000
*****
PH CALIB DATA
CALIB TIME:                        08:12:00
                                2008/06/01
OPERATOR NO:                       000

                                POINT 1
pH:                                4.003pH
POTENTIAL:                         177.3mV
TEMP:                              25.0c
*****
                                POINT 2
pH:                                6.864pH
POTENTIAL:                         8.0mV
TEMP:                              25.0c
*****
                                POINT 3
pH:                                9.182pH
POTENTIAL:                        -129.1mV
TEMP:                              25.0c
*****
CALIB RESULT
SLOPE 1:                           100.00%
E0 1:                              -0.0mV
SLOPE 2:                           100.00%
E0 2:                              -0.0mV
=====
```

Figure 28 Printout format

Software Operation


To install the software

Agilent offers the optional Electrochemical Data Collecting Software (G4390A) for communication between the 3200P pH Meter and a PC.

To install the optional software, insert the DVD into your PC drive, then follow the instructions displayed on the computer for installing. After installation, open the software from a desktop icon or from the Start menu.

To run software

- 1 Turn on the meter.
- 2 Connect the meter to the computer using the USB cable.
- 3 Start the software. The software will recognize meter model and type. and provide corresponding functions. Refer to the software documentation.

3 Troubleshooting and Maintenance

Troubleshooting Common Problems	50
General Troubleshooting Procedure	52
Meter Self-Diagnostics and Messages	56
Maintenance	59
Electrode cleaning	59
Electrode storage	59
Consumables and replacement parts	60

This chapter describes how to verify whether the meter is working properly and how to maintain the meter.

Troubleshooting Common Problems

Table 3 lists common problems, possible causes, and how to resolve them.

Note that improper operation of the electrode will lead to abnormal readings. During measurement, soak the electrode tip completely in solution. The location of the electrode in the solution should have good mobility.

Table 3 3200P pH Meter troubleshooting

Number	Failure mode	Failure cause	Solutions
1	No display after meter starts up. The display is not lit up.	1 Power adaptor not installed correctly. 2 The power supply does not meet with the requirements. 3 The power adaptor is damaged. 4 The power socket has poor contact. 5 The LCD is damaged (After start up for a few minutes, the user can hear a buzzing sound after pressing [On/Off] key, but there is no display on LCD.)	1 Reconnect. 2 Use the required power supply. 3 Replace the power adaptor. 4 Ensure there is a good contact. 5 Contact Customer Service.
2	Display is dim.		1 Check the power adapter.
3	No buzzing when pressing key	1 You may have pressed invalid keys under the current setup. 2 The buzzer has been damaged.	1 The buzzer only sounds when you press a valid key.
4	No response when pressing key	1 You may have pressed invalid keys under the current setup. 2 The key has been damaged.	1 Press the valid key for operation.
5	Meter does not display 25 °C when ATC probe is not connected with meter.	When the meter is not connected with a ATC probe, the meter will automatically switch to manual temperature setting.	Manually input the correct solution temperature.

Table 3 3200P pH Meter troubleshooting (continued)

Number	Failure mode	Failure cause	Solutions
6	The reading is not stable after a long time.	<div><div>1</div>The electrode has been damaged or aged.</div> <div><div>2</div>There is strong electrical signal interference source nearby (electrical leaking, strong electro-magnetic field, etc.).</div>	<div><div>1</div>Replace the electrode.</div> <div><div>2</div>Remove the electrical signal interference source. Lift or move the beaker away from the source of interference. Shield the meter and beaker from the electromagnetic field. Use the ground line. For the ground wire, one terminal is connected with the meter and the other terminal is connected with interference source.</div>

General Troubleshooting Procedure

During the use of the meter, there are many factors that may affect the measurement, including the electrode, the standard solutions used to calibrate electrode, the sample solution, the temperature during measurement, incorrect operation, or use, and so forth. When the measured results are significantly different from what was expected, first determine whether the meter itself or factors other than the meter caused the error. Follow the suggestions described below to diagnose the problem. Since there are multiple measuring methods, you will need to diagnose the meter based on your application conditions.

Diagnose whether the errors are caused by meter hardware, electrode calibration data error, parameter errors, or other errors, including the sample solution, solution preparation, buffer solutions used for calibration, and so forth.

For factors other than the meter, compare measurement results between known solutions. Put the electrode in different standard solutions. Check the corresponding potential (mV) value and pH value. Based on the comparison, judge whether the deviation is resulted from the electrode, the solution, or something else.

Check the temperature measurement

Because temperature measurements are used for all measurements, always ensure temperature measuring is functioning well.

Connect the ATC temperature diagnostic tool (5185-8390) to the meter. Turn on the meter to enter into measurement state. With the ATC temperature diagnostic tool connected, the meter should display a temperature reading between 49.0 to 51.0 °C. If yes, the meter is correctly measuring temperature. If the meter displays a temperature reading significantly different from 50 °C, there is a problem with the meter. Contact Agilent service.

Check the potential (mV) value

Disconnect the pH electrode. Install the short circuit plug (G4383-40000, shipped with the meter) into the pH/pX socket. Turn on the meter. Enter the pH measurement state. At this moment, the meter should display a potential (mV) reading between -0.5 to 0.5 mV (also called mV Zero). If yes, this indicates that the meter is functioning properly. If the potential (mV) reading significantly deviates from zero, then use the Adjust mV Zero function to adjust mV Zero. See [“To adjust mV zero”](#) on page 33 for more information. After adjustment, the potential (mV) reading should be near to zero. If the potential value (mV) reading cannot be adjusted near to zero by the mV Zero adjustment, this indicates a problem with the meter. Contact Agilent service.

Check the pH value

Disconnect the ATC probe. Set up the meter as manual temperature input. See [“To manually set the solution temperature”](#) on page 27. Set the solution temperature to 25 °C. When the potential (mV) is zero, the meter should display a pH reading among 6.99 to 7.01 pH. If yes, this indicates the electrode calibration data of the meter is acceptable. If not, there may be a problem with the electrode calibration data. Check the slope of electrode again. Press **[View]**, select **Last pH Calib** and press **[Enter]**. Now, view the calibration data and slope of last calibration. If the calibration data is severely out of the expected range (the meter displays an error message, see [Table 4](#) on page 57), this indicates a problem in the last electrode calibration, which may lead to a large measurement error. Recalibrate the electrode. Pay attention to the Buffer Group selection during calibration.

NOTE

In the measurement state, the meter may display an error message simultaneously with measurement data.

Restore factory default settings

If other troubleshooting does not find the problem, you can restore the meter to the factory default settings. This clears all user calibrations and settings so that you can verify that an erroneous calibration or setting is not the cause of the problem. After restoring the factory default settings, retest the meter.

To restore factory defaults (Set Default)

Select **Set Default** from the initial state display to restore all parameters to the factory defaults settings. The factory defaults are listed below.

- Clear mV Zero:

Any user-set mV zero value is cleared. With the short-circuit plug connected, check the mV reading. If the meter displays a non-zero mV value, adjust the mV zero. Failure to check and adjust the mV zero can lead to incorrect measurements.

- Modify Buffer Group data:

The Buffer Group of NIST standards includes three points, 4.01 pH, 7.00 pH, and 10.01 pH. The Buffer Group of DIN standards includes five points: 1.68 pH, 4.01 pH, 6.86 pH, 9.18 pH, and 12.45 pH. The Buffer Group of GB standards includes five points: 1.68 pH, 4.01 pH, 6.86 pH, 9.18 pH, and 12.45 pH. For example, the current Buffer Group is set up as GB standards.

- Modify pH Calibration data:

Three standard buffer solutions and the data are as listed below. By default, the electrode slope is set to 100% and E_0 is set to 0.0 mV.

- pH buffer 1: 4.003 pH/177.299 mV/25.0 °C
- pH buffer 2: 6.864 pH/ 8.046 mV/ 25.0 °C
- pH buffer 3: 9.182 pH/- 129.085 mV/ 25.0 °C
- Set manual temperature: 25 °C
- Set measuring method: **Continuous Mode**

- Set interval of **Timed Mode**: 10 minutes.
- Disable calibration message interval (**Calib Interval**)
- Set the meter to overwrite the oldest data with new data if the data storage becomes full.
- Set **Operator No.:** 000
- Set **Auto-Lock** condition in auto-lock reading mode: 0.1 pH
- Set **Auto-Lock Time** in auto-lock reading mode: 5 s
- Set pH resolution: 0.001 pH
- Set **Auto Power Off**: 0 (disabled)

Meter Self-Diagnostics and Messages

The meter supports self-diagnosis, which can find some common errors caused by meter, electrode, solution, or operation. It deals with these errors differently according to their severity.

There are two levels of severity: **Severe system errors** hinder further use of the meter or make the meter unable to finish necessary tasks. In this case, replacement or maintenance of the meter is required.

Minor errors, which are prompts or warnings, are caused by various reasons and can be corrected in several ways. However, you must pay attention to these errors to ensure measurement integrity and reliability. If you ignore these error messages the operation can continue.

Severe errors

When severe errors occur in the meter, the meter will shut down and display an error message similar to [Figure 29](#). In this case, the only operation that can be performed is to turn off the meter.


Figure 29 Severe error message

Error warning

Error warnings indicate improper electrode installation, solution problems, and electrode setup during use. Read these error warnings carefully. Proper maintenance and operation will reduce errors.

Figure 30 shows an example error warning that temperature is out of range.


Figure 30 Temperature warning

 is the error warning icon. The icon includes the error code. See Table 4 for details.

Table 4 Error warning codes

Code	Description	Solution	Memo
00	Conductivity measuring module error	Contact Customer Service.	Serious error
01	DO measuring module error	Contact Customer Service.	Serious error
02	Temperature measuring module error	Contact Customer Service.	Serious error
03	Data storage error	Contact Customer Service.	Serious error
20	Potential is out of range	Replace electrode.	–1999.9 to 1999.9 mV
21	pH/pX is out of range	Replace electrode.	–3.000 to 21.000 pH
22	Conductivity is out of range	Replace probe and solution.	0 to 2000 mS/cm
23	Resistivity is out of range	Replace probe and solution.	0 to 100 MΩ•cm

Table 4 Error warning codes (continued)

Code	Description	Solution	Memo
24	Temperature is out of range	Replace electrode and reduce solution temperature.	–6.0 to 120.0 °C
25	DO electric current is out of range	Replace electrode.	0 to 4000 nA
26	pH electrode slope is out of range	Replace electrode and recalibrate it.	80 to 120 %
27	Failed to recognize pH buffers	Replace electrode, setup proper buffer group and replace buffers.	
28	Temperature of pH buffer is out of range	Cool or heat buffer.	
29	Calibrate the same buffer repeatedly	Replace the buffer.	Due to incorrect operation
30	Number of pH buffers exceeds maximum.	Remove one or more buffers from the group.	5 buffers at most
31	pH buffers conflict with each other	Remove one or more buffers with neighboring pH values.	
32	The data storage is full.	Cover the previous data and store new data.	Do not delete all data
33	The internal clock has a low battery.	Set time manually.	
34	Number of customer-defined ions exceeds maximum	Delete one or more unnecessary customer-defined ions.	
35	The maximum ion mode number stored	Delete all storage data of certain ion mode.	

Maintenance

This section describes how to maintain and store the electrode.

Electrode cleaning

General electrode cleaning includes inorganic cleaning, organic cleaning, grease cleaning, protein precipitation cleaning and glass sensitive membrane regeneration. The type of cleaning needed depends on the contaminants and the electrode. After one or more cleaning procedures, rinse the outside of the electrode with distilled water. Siphon the reference filling solution away and add fresh solution. Repeat 2 or 3 times. Store the PH electrode in reference filling solution for at least one hour after cleaning. Store the pH electrode in the storage solution.

Electrode storage

When the electrode will not be used for 3 hours or more, store it in reference filling solution. For refillable electrodes, plug the filling hole. When the electrode is not used for a long time, put the electrode back in the box and store it a dry place at ambient temperature.

Consumables and replacement parts

Table 5 lists the consumables and replacement parts for the 3200P pH Meter.

Table 5 3200P pH Meter Consumables


Order No.	Model and name	Description
G4388A	3200SA Stirrer	Electrode holder and magnetic stirrer should be combined to stir solution with stable and precise speed with a large adjustable range.
G4389A	3200EA Electrode holder	It is used to support different probes and electrodes.
G4390A	Software package	3200P pH Meter can communicate with computer after this software package is installed (including software and communication cable).
5185-8389	Power adaptor	100–240 VAC, 1 A, 50/60 Hz
G4388-27000	Stirring bar	It is used with stirrer.
5190-3988	P3211 pH Combination electrode	Combination electrode, glass body and refillable
5190-3989	P3212 pH Combination electrode	Combination electrode, plastic body and non-refillable
5190-3990	P3311 pH Triode combination electrode	Triode combination electrode, with built-in ATC probe, glass body and refillable
5190-3991	P3111 pH Electrode	Glass body, single electrode
5190-3992	P3213 pH Combination electrode	Combination electrode, plastic body, flat sensitive glass membrane, refillable
5190-3993	P3214 pH Combination electrode	Combination electrode, plastic body, spear-tip sensitive glass, and non-refillable
5190-4003	R8111 Reference electrode	Glass body, ceramic junction, Ag/AgCl and single-junction
5190-3999	ORP8211 ORP electrode	Combination electrode, glass body and refillable
5190-3998	T7111 ATC Probe	Stainless Steel Temperature measurement range: 0–100 °C
5185-8390	ATC temperature diagnostic tool	Temperature troubleshooting
5185-8391	Conductivity diagnostic tool	Conductivity troubleshooting
G4383-40000	Short circuit plug	Used for zero potential calibration and troubleshooting.

Table 5 3200P pH Meter Consumables (continued)

Order No.	Model and name	Description
5190-0533	pH Buffer Kit (NIST Standard)	4.01, 7.00, 10.01 pH; Bottled Solution, 250 mL/bottle
5190-0534	pH Buffer Kit (GB Standard)	4.00, 6.86, 9.18 pH Bottled Solution, 3 X 250 mL/bottle
5190-0538	pH 4.00 Buffer (GB Standard)	Bottled Solution, 3 X 250 mL/bottle
5190-0535	pH 4.01 Buffer (NIST Standard)	Bottled Solution, 3 X 250 mL/bottle
5190-0539	pH 6.86 Buffer (GB Standard)	Bottled Solution, 3 X 250 mL/bottle
5190-0537	pH 7.00 Buffer (NIST Standard)	Bottled Solution, 3 X 250 mL/bottle
5190-0540	pH 9.18 Buffer (GB Standard)	Bottled Solution, 3 X 250 mL/bottle
5190-0536	pH 10.01 Buffer (NIST Standard)	Bottled Solution, 3 X 250 mL/bottle
5190-0545	Reference Filling Solution	Bottled Solution, 3 X 30 mL/bottle

Visit the Agilent website at <http://www.agilent.com/chem> for the latest consumables and replacement parts.

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